



Looking for
Dreve plastics in
the **USA** or **Canada**?

Available at
orthodontics.dreve-america.com

✉ sales@dreve-america.com

☎ +1 612-800-0021

Thermoforming blank made of ethyl vinyl acetate. Soft elastic with high recovery from deformation with a final hardness of 82 Shore A. Can be laminated and is ideally suited for sports mouthguards, positioners and duplication forms.

Please note:

- Drufosoft® transparent and colour are available in extra sizes according to inquiry
- Drufosoft® bi- and tricolour as well as colour mix are only available with diameter 120 mm

mm	colour	pieces	item No.
1	transparent	20	D4245
		50	D42451
1.5	transparent	20	D4252
		50	D42531
2	transparent	15	D4250
		50	D42511
3	transparent	10	D4248
		30	D42481
	neon-red	3	D42483X2
		10	D4248X2
	neon-green	3	D42483X3
		10	D4248X3
	neon-yellow	3	D42483X4
		10	D4248X4
	blue	3	D42483X5
		10	D4248X5
	yellow	3	D42483X6
		10	D4248X6
black	3	D42483X7	
	10	D4248X7	

mm	colour	pieces	item No.
3	white	3	D42483X8
		10	D4248X8
	pink	3	D42483X9
		10	D4248X9
	red	3	D42483X17
		10	D4248X17
	turquoise	3	D42483X18
		10	D4248X18
	gold	3	D42483X19
		10	D4248X19
	silver	3	D42483X20
		10	D4248X20
	blue-opak	3	D42483X21
		10	D4248X21
	green	3	D42483X22
		10	D4248X22
	purple	3	D42483X23
		10	D4248X23
	light blue	3	D42483X24
		10	D4248X24
bordeaux	3	D42483X25	
	10	D4248X25	
orange	3	D42483X26	
	10	D4248X26	
3	bicolour (yellow/black/yellow)	3	D4270X11
	bicolour (white/blue/white)		D4270X12
	bicolour (red/blue/red)		D4270X13
	tricolour (black/red/yellow)		D4275X14
	tricolour (blue/white/red)		D4275X15
	tricolour (red/green/blue)		D4275X16
	colour mix*		25
4	transparent	15	D42472
		75	D424901
5	transparent	15	D4225

* 1 blank of each colour and colour combination



Drufosoft® tricolour



Drufosoft® bicolour



Drufosoft® colour

Kombiplast

Highly transparent two layer thermoforming blank made of PETG and EVA. High wearing comfort occlusion splint through hard surface and soft inner side. Bonds with acrylate for the individual adjustment.

mm	pieces	item No.
2 (soft 1 : hard 1)	10	D420016
3 (soft 1.5 : hard 1.5)	10	D420017
4 (soft 1 : hard 3)	10	D420018



Occlusal splint made of Kombiplast

Drufosoft® pro

Transparent thermoforming blank made of ethyl vinyl acetate. Hard-flexible occlusion splints with high recovery from deformation with a final hardness of 95 Shore A. Can be laminated and is perfectly suitable for occlusion splints, positioners and bruxism splints.

mm	pieces	item No.
3	10	D4230



Occlusal splint made of Drufosoft® pro

Biolon

The thermoforming blank made of PETG (polyethylene terephthalate glycol) is supplied with an insulating film for the purpose of high transparency. Since it is tough-elastic and can be built up with acrylate, it is particularly suitable for aligner, bite and occlusion splints, as well as for orthodontic appliances.

mm	pieces	item No.
0.5	25	D420001
	50	D420002
0.625	20	D420019
	40	D420020
0.75	20	D420003
	40	D420004
1	25	D420005
	50	D420006
1.5	30	D420007
	60	D420008
2	20	D420009
	60	D420010
3	10	D420011
	50	D420012



Aligner made of Biolon