

# FotoDent® Materialmatrix



Drucker / Printer	Modell	denture	splint	IBT	guide	tray2	cast	model2 beige	model2 black	model2 white	gingiva	biobased model
Phrozen	Sonic XL 4K (2022) Qualified by Dreve	●	●	●	●	●	●	●	●	●	●	●
Ackretta	SOL	○		○	●	○	○	●	○	●	○	○
ASIGA	MAX UV	●		●	●	●	●	●	●	●	●	●
	PRO 4K	●		●	●	●	●	●	●	●	●	●
atum 3D	DLP-STATION 5-405											●
Carbon	M1			●			●				●	
	M2			●		●	●				●	
	M3			●		●	●				●	
	M3 Max			●		●	●				●	
Dekema	trix print2							●		●		
LuxCreo	iLux Pro Dental	○		○	○	○					○	○
Microlay	EVE	○										●
	VERSUS	○										●
MiiCraft	125	●			●		●	●	●	●	●	
Photocentric	Liquid Crystal OPUS	○		○	○							○
rapidshape	D series	●		●	●	●	●	●	●	●	●	●
	DII series	●		●	●	●	●	●	●	●	●	●
	D20+/D30+	●		●	●	●	●	●	●	●	●	●
	D90II	●										
Shining3D	AccuFab-L4D	○		○	○	○	○	○	○	○	○	○
Stratasys	Origin One	●										
W2P	SolFlex 350											●
	SolFlex 650											●

● Druckparameter erhältlich / Print parameter available

○ Druckparameter demnächst erhältlich / Print parameter available soon